

# HybriMore™ Hybridoma Culture Supplement

**HB01-1L**

V4.1

Store at 2-8 °C  
For Research Use Only

## ■ Introduction

Hybridoma cells, which are used for production of monoclonal antibodies, result from fusion of an antibody producing B-cell with a tumor cell, typically myeloma cells. Growth factors and serum are generally used to support the hybridoma development and achieve an optimal cell density and cloning efficiency.

**HybriMore™ Hybridoma Culture Supplement** is the first chemically defined growth promoting supplement especially for hybridoma cells. It may substantially increase the cloning efficiency of the newly fused hybridoma cells. It also significantly increases the successful rate of monoclonalization. The addition of **HybriMore™ Hybridoma Culture Supplement** shows no negative effect on the yield of secreted antibody.

## ■ Product Components

### **HybriMore™ Hybridoma Culture Supplement (HB01-1L)**

Hybridoma Culture Supplement  
User's manual

1 L dilute

1 bottle

## ■ Safety Information

Please wear gloves, lab coat and goggles while operating. Prevent contact product directly. In case of contacting, wash with large amount of water.

## ■ Storage

**HybriMore™ Hybridoma Culture Supplement** should be stored at 2-8 °C and shielded from light. Expiration date is labeled on the bottle or box.

## ■ Instruction

1. Reconstitute Hybridoma Culture Supplement by 1 mL ddH<sub>2</sub>O.
2. Add the reconstituted Hybridoma Culture Supplement into 1L regular culture media, such as DMEM or RPMI-1640 based media with 10% fetal bovine serum and antibiotics.
3. Sterile the reconstituted Hybridoma Culture Supplement by filtering through 0.22 µm filters. The reconstituted Hybridoma Culture Supplement should be used up within 1 month.
4. The Hybridoma Culture Supplement containing media can be stored as the regular hybridoma culturing media. The addition of Hybridoma Culture Supplement will not alter the shelf life of prepared media.
5. The prepared media can be directly used for culturing the newly fused hybridoma cells, or the monoclonization experiment.
6. The recommend situations of Hybridoma Culture Supplement addition:
  - (1) Hybridoma Culture Supplement can be added to HAT-media during the initial selection of hybridoma cells.
  - (2) Hybridoma Culture Supplement can be added to cell culture media for the monoclonization of hybridoma cells.

## ■ Related Visual Protein Products

CellCount™ MTT Assay Kit	CC01-11	1,000 tests
CellCount™ MTT Assay Kit (with MTT Solvent)	CC01-12	1,000 tests
CellCount™ MTT Assay Kit	CC01-51	5,000 tests
CellCount™ MTT Assay Kit (with MTT Solvent)	CC01-52	5,000 tests
CellCount™ Cell Counting Kit-8 (CCK-8)	CC02-01	500 tests
CellCount™ Cell Counting Kit-8 (CCK-8)	CC02-06	3,000 tests
CellCount™ Cell Counting Kit-8 (CCK-8)	CC02-20	10,000 tests
CellCount™ Cell Counting Kit-8 (CCK-8 powder)	CC03-01	500 tests
CellCount™ Cell Counting Kit-8 (CCK-8 powder)	CC03-10	5,000 tests
CytoMore™ Cell Rescue Supplement	CT01-1BT	1 bottle
NeuronMore™ Neural Cell Culture Supplement	CT02-1L	1 bottle
ImmunoFast™ Adjuvant	IF01-4N	4 reactions
ImmunoFast™ Adjuvant	IF01-20N	20 reactions
Trypan Blue Solution (0.4%)	TPB01-100	100 mL