

CellCount™ Cell Counting Kit-8 (CCK-8)

CC02-01/CC02-06/CC02-20

V3.2

Store at -20 °C
For Research Use Only

■ Introduction

The **CellCount™ Cell Counting Kit-8 (CCK-8)** is a non-radioactive assay kit that uses a water-soluble tetrazolium compound to measure cell viability. When the compound is added directly to cells, an electron carrier converts it into a water-soluble orange formazan dye, the amount of which is directly proportional to the number of living cells. The **CellCount™ Cell Counting Kit-8 (CCK-8)** offers a fast and convenient, ready-to-use method for cell proliferation and cytotoxicity studies. It has excellent stability and little cytotoxicity for living cells during long-term observation. The detection sensitivity of **CellCount™ Cell Counting Kit-8 (CCK-8)** is significantly higher than other similar assays using tetrazolium salts such as MTT, XTT, MTS or WST-1.

■ Product Components

| | | | |
|---|---------|------|---------------------|
| CellCount™ Cell Counting Kit-8 (CCK-8) (CC02-01) | | | 500 tests |
| Cell Counting Kit-8 (CCK-8) | CC02-01 | 5 mL | 1 bottle |
| User's manual | | | |
| CellCount™ Cell Counting Kit-8 (CCK-8) (CC02-06) | | | 3,000 tests |
| Cell Counting Kit-8 (CCK-8) | CC02-01 | 5 mL | 6 bottles |
| User's manual | | | |
| CellCount™ Cell Counting Kit-8 (CCK-8) (CC02-20) | | | 10,000 tests |
| Cell Counting Kit-8 (CCK-8) | CC02-01 | 5 mL | 20 bottles |
| User's manual | | | |

■ Safety Information

Please use personal protective equipment such as gloves, lab coat, and goggles when handling. Avoid direct contact with the product content. In the event of contact, wash the affected area with a large amount of water.

■ Storage

Store **CellCount™ Cell Counting Kit-8 (CCK-8)** at -20 °C and shielded from light. Expiration date is labeled on the bottle or box. Aliquot upon using and store working solutions at 2-8 °C for up to 3 months. Avoid repeated freeze-thaw cycle.

■ Materials needed but not provided

1. 96 well plate with clear bottom
2. Humidified incubator (e.g. CO₂ incubator)
3. Plate Reader capable of measuring absorbance in the region of 450 nm
4. 1% w/v SDS (Stop solution)

■ Instruction

NOTE: Sterilize CCK-8 solution by membrane filtration if necessary.

A. Cell Number Determination

1. Add 100 µL cell suspension containing known numbers of viable cells into a 96 well plate for a calibration curve.
2. Add 100 µL of the cell suspension to other wells then culture the cells in humidified incubator at 37 °C for 24 hours.
3. Thaw the CCK-8 solution on the benchtop before use.

NOTE: Thawed CCK-8 solution can be stored at 2-8 °C for 3 months.

NOTE: Repeated freeze-thaw will result in increased background signal affecting assay sensitivity.

4. Add 10 µL of the CCK-8 solution to each well. Avoid introducing bubbles to minimize interference to optical density reading.
5. Incubate the plate in humidified incubator at 37 °C for 1-4 hours.

NOTE: The incubation time can range from less than 1 to more than 4 hours, depending on cell type and cell numbers. Optimize the incubation time for each experiment.

6. Measure the absorbance at or near 450 nm on a plate reader. The optical density will be stable for 2 days by adding 10 µL of 1% w/v SDS (Stop solution) to each well.

B. Cell Proliferation Assay and Cytotoxicity Assay

1. Prepare a cell suspension (10,000-500,000 cells/ml, depending on application and cell type) using an appropriate culture medium.
2. Add 100 μ L of the cell suspension to each well of a 96-well plate then culture the cells in humidified incubator at 37 °C for 24 hours.
3. Add test articles at a volume of \sim 10 μ L into the test well.
4. Incubate the plate for an appropriate length of time (e.g. 6, 12, 24 or 48 hours) in the humidified incubator.
5. Thaw the CCK-8 solution on the benchtop before use.

NOTE: Thawed CCK-8 solution can be stored at 2-8 °C for 3 months.

NOTE: Repeated freeze-thaw will result in increased background signal affecting assay sensitivity.

6. Add 10 μ L of the CCK-8 solution to each well. Avoid introducing bubbles to minimize interference to optical density reading.
7. Incubate the plate in humidified incubator at 37 °C for 1-4 hours.

NOTE: The incubation time can range from less than 1 to more than 4 hours, depending on cell type and cell numbers. Optimize the incubation time for each experiment.

8. Measure the absorbance at or near 450 nm on a plate reader. The optical density will be stable for 2 days by adding 10 μ L of 1% w/v SDS (Stop solution) to each well.

Troubleshooting

| Problem | Possible cause | Remedy |
|--|---|--|
| The optical density reading is higher than expected | High background | Avoid repeated thawing and freezing the CCK-8 solution |
| | Too many cells in a well | Decrease the cell numbers per well or decrease the incubation time of CCK-8 solution |
| The optical density reading is lower than expected or remains the same as the blank well | Very few cells in a well | Increase the cell numbers per well or increase the incubation time of CCK-8 solution |
| | CCK-8 solution is not added to the well | Add the CCK-8 solution to each well |
| | The substances cause cell death | Decrease the concentration of substances |

■ Related Visual Protein Products

| | | |
|---|-----------|--------------|
| CellCount™ MTT Assay Kit | CC01-11 | 1,000 tests |
| CellCount™ MTT Assay Kit (with MTT Solvent) | CC01-12 | 1,000 tests |
| CellCount™ MTT Assay Kit | CC01-51 | 5,000 tests |
| CellCount™ MTT Assay Kit (with MTT Solvent) | CC01-52 | 5,000 tests |
| CellCount™ Cell Counting Kit-8 (CCK-8 powder) | CC03-01 | 500 tests |
| CellCount™ Cell Counting Kit-8 (CCK-8 powder) | CC03-10 | 5,000 tests |
| CytoMore™ Cell Rescue Supplement | CT01-1BT | 1 bottle |
| NeuronMore™ Neural Cell Culture Supplement | CT02-1L | 1 bottle |
| HybriMore™ Hybridoma Culture Supplement | HB01-1L | 1 bottle |
| ImmunoFast™ Adjuvant | IF01-4N | 4 reactions |
| ImmunoFast™ Adjuvant | IF01-20N | 20 reactions |
| Trypan Blue Solution (0.4%) | TPB01-100 | 100 mL |